

Estimating Milling Time

Milling rate depends on:

- Beam Current
- Volume to be milled
- A material-unique constant (volume/dose, 'V/d'):

$$t_{\text{mill}}[\text{min}] = \frac{(1\text{min}/60\text{s}) * V[\mu\text{m}^3]}{(V/d) [\mu\text{m}^3/\text{nC}] * I[\text{nA}]} \approx \frac{(\text{depth}[\mu\text{m}])^3}{73 * (V/d) * I[\text{nA}]}$$


Table 5-7 Material Sputter Rates at 30 kV

Material	Volume per Dose [$\mu\text{m}^3 / \text{nC}$]	Material	Volume per Dose [$\mu\text{m}^3 / \text{nC}$]
C	0.18	Au	1.50
Si	0.27	MgO	0.15
Al	0.30	SiO ₂	0.24
Ti	0.37	Al ₂ O ₃	0.08
Cr	0.10	TiO	0.15
Fe	0.29	Si ₃ N ₄	0.20
Ni	0.14	TiN	0.15
Cu	0.25	Fe ₂ O ₃	0.25
Mo	0.12	GaAs	0.61
Ta	0.32	Pt	0.23
W	0.12	PMMA	0.40

Volume rules of thumb

Cross-section milling generates a wedge (min. volume to expose face to e⁻ beam): $V_{\text{mill}} = 1/2 l * w * d$

- $w \approx l$
- $l \geq d * \tan(52)$
- $V_{\text{mill}} \approx 0.82 d^3$



Example calculations for various materials and depths

Z	(V/d)	q [$\mu\text{m}^3/\text{s}$]	t_{mill} [min]	time (d=5 μm)	time t(d=20 μm)	time (d=100 μm)
Au	1.50	97.5	$d^3/7135$	1 s	1 min	2.5 hr
Al	0.30	19.5	$d^3/1427$	5 s	5.5 min	12 hr
Al ₂ O ₃	0.08	5.2	$d^3/381$	20 s	21 min	44 hr

Estimating Milling Time

The listed sputter rates have a roughly ~linear dependence on the hardness. So, for materials not listed, but which have known hardnesses, the sputtering time can be estimated by:

$$t_{\text{mill}}[\text{min}] \approx \frac{(\text{depth}[\mu\text{m}])^3}{4.83 * (10 - \text{Mohs}) * I[\text{nA}]}$$




Table 5-7 Material Sputter Rates at 30 kV


Material	Volume per Dose [μm ³ / nC]	Material	Volume per Dose [μm ³ / nC]
C	0.18	Au	1.50
Si	0.27	MgO	0.15
Al	0.30	SiO₂	0.24
Ti	0.37	Al₂O₃	0.08
Cr	0.10	TiO	0.15
Fe	0.29	Si₃N₄	0.20
Ni	0.14	TiN	0.15
Cu	0.25	Fe₂O₃	0.25
Mo	0.12	GaAs	0.61
Ta	0.32	Pt	0.23
W	0.12	PMMA	0.40

Estimating Milling Time

For unknown materials (or for a more precise value), a short preliminary test can determine milling rate:



- Mill a small square ($\sim 5 \times 5 \mu\text{m}$) in a test region for 60s
- Measure the depth achieved
- $(V/d)[\mu\text{m}^3/\text{nA}/\text{min}] = 25 * \text{depth}[\mu\text{m}]/I[\text{nA}]$

- *In the example to the right, using the Moh's value for Ag (2.5-3) yields a very similar milling rate compared to what we see experimentally (~ 31 vs $\sim 35 \mu\text{m}^3/\text{min}/\text{nA}$)*
- *Conversely, the predicted rate for SiC (Moh's 9-9.5) is significantly lower than what we actually see (~ 3 vs $10 \mu\text{m}^3/\text{min}/\text{nA}$)*



Oblique Milling (ie milling w/o tilting the sample perpendicular to the ion beam). Useful for:

- Samples that aren't easily tilted
- Analysis of very fine films – angled cut elongates the exposed cross-section



Requires measuring the thickness of an exposed film via this method requires **manual** conversion. Using the length acquired without any tilt correction, use the following:

$$l_{act} = \frac{l_{meas}}{\sin(52^\circ)} = 1.27 * l_{meas} \quad l_{film} = \frac{l_{meas}}{\tan(52^\circ)} = 0.78 * l_{meas}$$